
Sea sand, Hi-AR[®]

GRM3061

Product Identifier

CAS No.	:	7631-86-9
EC No.	:	231-545-4
Molecular Formula	:	SiO ₂
Molecular Weight	:	60.08
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to light brown to grey solid or powder
Solubility	:	Insoluble in water.
Particle size	:	85.00 - 100.00% passing through 20 mesh (~850 µm)
Loss on drying (at 105°C, 2 hr)	:	<= 0.20%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg