

**GRM3104 Ammonium phosphomolybdate**

**Product Number Packing**

---

GRM3104 : 25G

**Product Information**

---

Product Code : GRM3104  
Product Name : Ammonium phosphomolybdate  
Molecular Formula :  $H_{12}Mo_{12}N_3O_{40}P$   
Molecular Weight : 1876.35  
CAS No. : 12026-66-3  
EC No. : 234-714-0  
Shelf Life : 4 years

**Technical Specification**

---

Appearance : Yellow crystals or powder  
Solubility : 33.3 mg soluble in 1 mL of alkali hydroxide  
Assay (Lead Nitrate Titration) : 98.00 - 102.00 %

**Risk and Safety Information**

---

Storage Temperature(°C) : Store below 30°C

**Transport Information**

---

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods