
Bovine serum albumin pH 7.0, Microbiological grade

GRM3115

Product Identifier

CAS No.	:	9048-46-8
EC No.	:	232-936-2
HS Code	:	3502 90 00
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to light brown crystals or powder or flakes prepared by cold ethanol plasma fractionation
Solubility	:	33.3 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	6.5 - 7.5
Loss on drying (at 105°C, 2 hr)	:	<= 5.0%
Assay (Protein, on dry basis)	:	min. 95.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1

Other Information

Prepared by cold ethanol plasma fractionation. An exceptional material for use in tissue culture work.