
Bovine serum albumin, Cohn Fraction V powder

GRM3151

Product Identifier

| | | |
|------------|---|------------|
| CAS No. | : | 9048-46-8 |
| EC No. | : | 232-936-2 |
| HS Code | : | 3502 90 00 |
| Storage | : | 2 - 8°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|---------------------------------|---|---|
| Appearance | : | White to light brown crystals or powder prepared by cold ethanol plasma fractionation |
| Solubility | : | 33.3 mg soluble in 1 mL of water |
| pH (1% in water at 25°C) | : | 6.5 - 7.5 |
| Loss on drying (at 105°C, 2 hr) | : | <= 5.0% |
| Assay (Protein, on dry basis) | : | min. 96.00% |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | Not available |
| WGK | : | 1 |

Other Information

Prepared by cold ethanol plasma fractionation. Rich in lipids and cofactors