

---

## Azure A, Hi-Cert<sup>™</sup>

**GRM318**

### Product Identifier

---

CAS No.	:	531-53-3
EC No.	:	208-510-7
C.I. No.	:	52005
Molecular Formula	:	C <sub>14</sub> H <sub>14</sub> N <sub>3</sub> SCl
Molecular Weight	:	291.80
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Green to green with brown cast to black powder
Solubility	:	20 ppm in water yields clear solution
Absorption maxima	:	620 - 634 nm (20 ppm in water)
Loss on drying (at 110°C, 2 hr)	:	<= 8.0%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	SP5660000
WGK	:	3