

## Canada balsam, Synthetic

**GRM322**


### Product Identifier

CAS No.	:	8007-47-4
EC No.	:	232-362-2
Synonym	:	Balsam canada, Synthetic
HS Code	:	3301 30 99
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Colourless to yellow to greenish-yellow viscous liquid
Solubility	:	1 mL miscible in 1 mL of chloroform
Refractive index (n <sub>20/D</sub> )	:	1.520 - 1.540
Density (at 25°C)	:	0.910 - 0.990 g/mL

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H226
Precautionary Statement(s)	:	P210- P241- P280- P303+P361+P353- P370 + P378- P501
UN No.	:	1993
Class	:	3
Packing Group	:	III
RTECS	:	CP2352500
WGK	:	2