
Cholic acid sodium salt, Hi-LR™

GRM3269

Product Identifier

CAS No.	:	361-09-1
EC No.	:	206-643-5
Molecular Formula	:	C ₂₄ H ₃₉ NaO ₅
Molecular Weight	:	430.60
HS Code	:	2918 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white solid or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Water (K.F.)	:	<= 5.00%
Assay (GC/NT, on anhydrous basis)	:	>= 98.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1