

GRM3277 Cinnamaldehyde**Product Number Packing**

GRM3277 : 500ML

Product Information

Product Code : GRM3277
Product Name : Cinnamaldehyde
Synonym : 3-Phenylprop-2-enal
Molecular Formula : C₉H₈O
Molecular Weight : 132.16
CAS No. : 104-55-2
EC No. : 203-213-9
HS Code : 2912 29 00
Shelf Life : 4 years

Technical Specification

Appearance : Colourless to yellow liquid
Solubility : 1 mL miscible in 1 mL of alcohol
FTIR (Liquid film) : Matches with the standard pattern
Refractive index (n_{20/D}) : 1.614 - 1.623
Density (at 25°C) : 1.045 - 1.055 g/mL
Assay (GC) : min. 98.00 %

GHS Safety Information

Hazard Statement(s) : H315-H317-H319-H335
Precautionary Statement(s) : P261-P280-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37/38-43
S-Phrase(s) : 26-36/37
WGK : 1
RTECS : GD6475000
Flash Point(°F) : 159.8 °F
Flash Point(°C) : 71 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**



Product Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein