

---

## Carboxymethylcellulose sodium salt, High viscosity (CMC)

**GRM329**

---

### Product Identifier

CAS No.	:	9004-32-4
HS Code	:	3912 31 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to yellow or light beige powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	6.50 - 8.50
FTIR	:	Matches with the standard pattern
Loss on drying (at 105°C, 2 hr)	:	<= 10.00%
Viscosity (1% in water at 25°C)	:	1500 - 3000 cP

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1