
Cesium chloride, Hi-AR™

GRM330

Product Identifier

CAS No.	:	7647-17-8
EC No.	:	231-600-2
Molecular Formula	:	CsCl
Molecular Weight	:	168.36
HS Code	:	2827 39 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white deliquescent crystals or powder or beads or chunks
Solubility	:	1000 mg soluble in 1 mL of water
Copper (Cu)	:	<= 0.003%
Iron (Fe)	:	<= 0.0003%
Lead (Pb)	:	<= 0.0001%
Magnesium (Mg)	:	<= 0.0005%
Sulfate (SO ₄)	:	<= 0.002%
Assay (AT/EDTA Titration)	:	99.50 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	FK9625000
WGK	:	1