

Cobalt metal powder

GRM3313


Product Identifier

CAS No.	:	7440-48-4
EC No.	:	231-158-0
EC Index No	:	027-001-00-9
Atomic Formula	:	Co
Atomic Weight	:	58.93
HS Code	:	8105 20 20
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Orange to dark grey to black powder
Solubility	:	33.3 mg soluble in 1 mL of dilute nitric acid
Assay (EDTA Titration)	:	min. 99.50%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H317- H334- H413
Precautionary Statement(s)	:	P261- P280- P342 + P311
UN No.	:	3089
Class	:	4.1
Packing Group	:	II
RTECS	:	GF8750000
WGK	:	NWG