
Cellulose microcrystalline, For column chromatography

GRM333

Product Identifier

CAS No.	:	9004-34-6
EC No.	:	232-674-9
HS Code	:	3912 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white powder
Solubility	:	Insoluble in water
pH (10% aqueous slurry at 25°C)	:	5.0 - 7.5
Identification test	:	Passes
Loss on drying (at 105°C, 2 hr)	:	<= 6.0%
Arsenic (As)	:	<= 0.002%
Heavy metals (as Pb)	:	<= 0.001%
Sulphated ash	:	<= 0.2%
Water soluble substances	:	<= 0.2%
Assay (FAS Titration, on dry basis)	:	97.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	FJ5691460
WGK	:	1