
Diatomaceous earth

GRM3380


Product Identifier

CAS No.	:	61790-53-2
HS Code	:	2512 00 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to tan to grey powder
Solubility	:	Insoluble in water
Loss on drying (at 105°C, 2 hr)	:	<= 10.00%
Loss on ignition (at 800°C)	:	<= 7.00%
Nonsiliceous substances	:	min. 25.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319- H335- H373
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3

Other Information

Highly insoluble inert material for use as a filtration aid.