
Clove oil

GRM340


Product Identifier

CAS No.	:	8000-34-8
HS Code	:	3301 29 21
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow to yellow-green to brown viscous liquid
Solubility	:	33.3 mg soluble in 1 mL of methanol
Refractive index (n 20/D)	:	1.5320 - 1.5345
Density (at 25°C)	:	1.035 - 1.060 g/mL
Assay (Eugenol)	:	min. 85.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	GF6900000
WGK	:	2