
EDTA disodium magnesium salt tetrahydrate

GRM3445

Product Identifier

CAS No.	:	29932-54-5
EC No.	:	238-372-3
Molecular Formula	:	$C_{10}H_{12}MgN_2Na_2O_8 \cdot 4H_2O$
Molecular Weight	:	430.56
Synonym	:	Ethylenediaminetetraacetic acid disodium magnesium salt tetrahydrate
HS Code	:	2922 49 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Water (K.F.)	:	10.00 - 20.00%
Assay (EDTA Titration)	:	99.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	2