
o-Cresolphthalein, Indicator grade

GRM345

Product Identifier

CAS No.	:	596-27-0
EC No.	:	209-881-8
Molecular Formula	:	C ₂₂ H ₁₈ O ₄
Molecular Weight	:	346.38
HS Code	:	2932 20 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to tan powder
Solubility	:	20 ppm in ethanol yields clear solution
Visual pH Transition	:	pH 8.0 (Colourless) - pH 9.8 (Pink/Violet Red)
Loss on drying (at 110°C, 2 hr)	:	<= 1.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	SM8390000
WGK	:	3