

## N,N'-Methylenebis(acrylamide), Hi-LR™

**GRM3489**


### Product Identifier

CAS No.	:	110-26-9
EC No.	:	203-750-9
Molecular Formula	:	C <sub>7</sub> H <sub>10</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight	:	154.17
Synonym	:	Bis-acrylamide
HS Code	:	2924 19 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to almost white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of methanol giving clear to hazy solution
FTIR	:	Matches with the standard pattern
Assay (GC/Nitrogen basis)	:	min. 98.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H332
Precautionary Statement(s)	:	P301+P312+P330
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AS3678000
WGK	:	2