
Lithium bromide

GRM3535


Product Identifier

CAS No.	:	7550-35-8
EC No.	:	231-439-8
Molecular Formula	:	LiBr
Molecular Weight	:	86.85
HS Code	:	2827 59 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white hygroscopic crystals or powder or granules or chunks
Solubility	:	100 mg soluble in 1 mL of water
Assay (AT)	:	99.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H317- H319
Precautionary Statement(s)	:	P261- P264.- P280- P301 + P312- P302 + P352- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	OJ5755000
WGK	:	1