

## 3-Mercaptopropionic acid

**GRM3550**


### Product Identifier

|                   |   |  |
|-------------------|---|--|
| CAS No.           | : | 107-96-0                                       |
| EC No.            | : | 203-537-0                                      |
| Molecular Formula | : | C <sub>3</sub> H <sub>6</sub> O <sub>2</sub> S |
| Molecular Weight  | : | 106.14   |
| HS Code           | : | 2930 90 99                                     |
| Storage           | : | Below 30°C                                     |
| Shelf life        | : | 4 years  |

### Technical Specification

|            |   |                                |
|------------|---|--------------------------------|
| Appearance | : | Colorless liquid.              |
| Solubility | : | 1 mL miscible in 1 mL of water |
| Assay (GC) | : | min. 99.00%                    |

### Safety Information

|                            |   |   |
|----------------------------|---|---|
| Hazard Pictogram(s)        | : |  |
| Signal Word                | : | Danger  |
| Hazard Statement(s)        | : | H301- H314  |
| Precautionary Statement(s) | : | P280- P301 + P310- P305+P351+P338- P310   |
| UN No.                     | : | 2922  |
| Class                      | : | 8   |
| Packing Group              | : | II  |
| RTECS                      | : | UF5270000   |
| WGK                        | : | 1   |