

Melamine

GRM3588


Product Identifier

CAS No.	:	108-78-1
EC No.	:	203-615-4
EC Index No	:	613-345-00-2
Molecular Formula	:	C ₃ H ₆ N ₆
Molecular Weight	:	126.12
Synonym	:	2,4,6-Triamino-1,3,5-triazine
HS Code	:	2933 61 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder
Solubility	:	10 mg soluble in 10 mL of water
FTIR	:	Matches with the standard pattern
Assay (NT/HPLC)	:	>= 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H351- H373
Precautionary Statement(s)	:	P201- P260- P280- P308 + P313- P314- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	OS0700000
WGK	:	2