

m-Nitrobenzenesulphonic acid sodium salt

GRM3638


Product Identifier

CAS No.	:	127-68-4
EC No.	:	204-857-3
EC Index No	:	609-048-00-2
Molecular Formula	:	C ₆ H ₄ NNaO ₅ S
Molecular Weight	:	225.15
HS Code	:	2904 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Off-white to yellow crystals or powder or chunks or flakes
Solubility	:	33.3 mg soluble in 1 mL of water (clear to slightly hazy)
FTIR	:	Matches with the standard pattern
Assay	:	95.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H317- H319
Precautionary Statement(s)	:	P280- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DB7195000
WGK	:	1