

**GRM365 Cadmium carbonate****Product Number Packing**

GRM365 : 500G

**Product Information**

Product Code : GRM365  
Product Name : Cadmium carbonate  
Molecular Formula :  $CdCO_3$   
Molecular Weight : 172.42  
CAS No. : 513-78-0  
EC No. : 208-168-9  
Shelf Life : 4 years

**Technical Specification**

Appearance : White powder or granules or chunks  
Solubility : 33.3 mg soluble in 1 mL of dilute acids  
Assay (Cd, EDTA Titration) : 62.00 - 67.40%

**GHS Safety Information**

Hazard Statement(s) : H302-H312-H332-H410  
Precautionary Statement(s) : P273-P280-P501  
Signal Word : Warning  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 20/21/22-50/53  
S-Phrase(s) : 60-61  
WGK : 3  
RTECS : FF9320000  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

UN No. : 3077  
Class : 9  
Packaging Group : 3  
Marine Pollutant : Yes  
ADR/RID : 3077 9/PG 3  
IMDG : 3077 9/PG 3  
IATA : 3077 9/PG 3