

Paraformaldehyde

GRM3660


Product Identifier

CAS No.	:	30525-89-4
Molecular Formula	:	HO(CH ₂ O) _n H
Molecular Weight	:	30.03 (as monomer)
HS Code	:	2912 60 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of 1M NaOH
FTIR	:	Matches with the standard pattern
Assay (H ₂ SO ₄ Titration)	:	95.00 - 106.50%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H228- H302+H332- H315- H317- H318- H335- H341- H350
Precautionary Statement(s)	:	P201- P210- P280- P305+P351+P338- P370 + P378- P403 + P233
UN No.	:	2213
Class	:	4.1
Packing Group	:	III
RTECS	:	RV0540000
WGK	:	2