
Prussian blue, Water soluble

GRM3665

Product Identifier

CAS No.	:	12240-15-2
Molecular Formula	:	$C_6Fe_2KN_6 \cdot xH_2O$
Molecular Weight	:	306.89 (anhydrous basis)
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Dark blue to very dark blue powder
Solubility	:	10ppm in water yields clear solution
Absorption maxima	:	253 - 263 nm (10 ppm in water)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1