
Silica gel GF 254, For TLC

GRM367

(contains CaSO₄ binder and fluorescent indicator)

Product Identifier

CAS No.	:	112926-00-8
EC No.	:	231-545-4
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white fine powder or crystals
Solubility	:	Insoluble in water
pH (10% aqueous suspension at 25°C)	:	6.50 - 7.50
Chloride (Cl)	:	<= 0.02%
Iron (Fe)	:	<= 0.02%
Suitability for TLC	:	Passes test
CaSO ₄ content (EDTA Titration)	:	12.50 - 13.50%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg