

## Stannous (II) sulphate

**GRM3726**


### Product Identifier

CAS No.	:	7488-55-3
EC No.	:	231-302-2
Molecular Formula	:	SnSO <sub>4</sub>
Molecular Weight	:	214.77
HS Code	:	2833 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to off-white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of dilute sulphuric acid
Assay (Iodimetry)	:	>= 96.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H290- H315- H317- H319- H332- H335- H341- H373- H410
Precautionary Statement(s)	:	P260- P273- P280- P312- P403 + P233- P501
UN No.	:	3260
Class	:	8
Packing Group	:	III
RTECS	:	Not available
WGK	:	3