

RM3731 Syringic acid**Product Number Packing**

RM3731 : 25G

Product Information

Product Code : RM3731
Product Name : Syringic acid
Synonym : 3,5-Dimethoxy-4-hydroxybenzoic acid, 4-Hydroxy-3,5-dimethoxy-benzoic acid
Molecular Formula : C₉ H₁₀ O₅
Molecular Weight : 198.17
CAS No. : 530-57-4
EC No. : 208-486-8
HS Code : 2918 99 90
Shelf Life : 4 years

Technical Specification

Appearance : White to greyish-beige to tan powder.
Solubility : Soluble in alcohol
Melting point : 205 - 209°C
FTIR (KBr disc) : Matches with the standard pattern.
Assay (NaOH titration) : min. 97.00 %

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305 + P351 + P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37/38
S-Phrase(s) : 26-36
WGK : 3
RTECS : DH2090000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods



Product Information

IATA

: Not Dangerous Goods

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein