
Fluorescein sodium salt

GRM374

Product Identifier

CAS No.	:	518-47-8
EC No.	:	208-253-0
C.I. No.	:	45350
Molecular Formula	:	C ₂₀ H ₁₀ Na ₂ O ₅
Molecular Weight	:	376.27
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	Light red to dark brown or greenish-brown crystals or powder or solid
Solubility	:	10 ppm in buffer pH 8 yields clear solution
Absorption maxima	:	488 - 491 nm (10 ppm in buffer pH 8)
	:	

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	LM5425000
WGK	:	1