

GRM3766 Triethyl phosphate**Product Number Packing**

GRM3766 : 500ML

Product Information

Product Code : GRM3766
Product Name : Triethyl phosphate
Synonym : Phosphoric acid triethyl ester; TEP
Molecular Formula : C₆H₁₅O₄P
Molecular Weight : 182.15
CAS No. : 78-40-0
EC No. : 201-114-5
HS Code : 2918 29 00
EC Index No. : 015-013-00-7
Shelf Life : 4 years

Technical Specification

Appearance : Colourless liquid
Solubility : 33.3 mg soluble in 1 mL of methanol
FTIR (Liquid film) : Matches with the standard pattern
Refractive index (n 20/D) : 1.3980 - 1.4080
Density (at 25°C) : 1.067 - 1.077 g/mL
Assay (GC) : min. 99.00%

GHS Safety Information

Hazard Statement(s) : H302
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22
S-Phrase(s) : 25
WGK : 1
RTECS : TC7900000
Flash Point(°F) : 266 °F
Flash Point(°C) : 130 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No



Product Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods