

## Tri-n-butyl phosphate

**GRM3771**


### Product Identifier

CAS No.	:	126-73-8
EC No.	:	204-800-2
EC Index No	:	015-014-00-2
Molecular Formula	:	C <sub>12</sub> H <sub>27</sub> O <sub>4</sub> P
Molecular Weight	:	266.31
HS Code	:	2919 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Colourless to light yellow liquid
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Refractive index (n 20/D)	:	1.4190 - 1.4290
Density (at 25°C)	:	0.974 - 0.984 g/mL
Assay (GC)	:	min. 99.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H351- H412
Precautionary Statement(s)	:	P273- P301+P312+P330
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TC7700000
WGK	:	2