

2-Chlorobenzoic acid, Hi-LRTM

GRM3985


Product Identifier

CAS No.	:	118-91-2
EC No.	:	204-285-4
Molecular Formula	:	C ₇ H ₅ ClO ₂
Molecular Weight	:	156.57
Synonym	:	o-Chlorobenzoic acid
HS Code	:	2916 39 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow or grey crystals or powder or chunks
Solubility	:	100 mg soluble in 1mL of alcohol
FTIR	:	Matches with the standard pattern
Melting range	:	137 - 144°C
Assay (NaOH Titration/GC)	:	97.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319
Precautionary Statement(s)	:	P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DG4976000
WGK	:	2