

**GRM4020****L-Amino acids kit****Product Number****Packing**

GRM4020	: 1KT
GRM4020	: 5KT

**Product Information**

Product Code	: GRM4020
Product Name	: L-Amino acids kit
HS Code	: 3822 00 00
Shelf Life	: 4 years

**Technical Specification**

Description	: The kit contains the following 22 individual amino acids, each having 1 gm lyophilized powder
L-Alanine	: Corresponding HiMedia Code : GRM036
L-Arginine monohydrochloride	: Corresponding HiMedia Code : GRM039
L-Asparagine monohydrate	: Corresponding HiMedia Code : GRM041
L-Aspartic acid	: Corresponding HiMedia Code : GRM043
L-Cysteine	: Corresponding HiMedia Code : GRM883
L-Cystine, Free base	: Corresponding HiMedia Code : GRM047
L-Glutamic acid	: Corresponding HiMedia Code : GRM048
L-Glutamine	: Corresponding HiMedia Code : RM049
Glycine	: Corresponding HiMedia Code : GRM199
L-Histidine monohydrochloride monohydrate	: Corresponding HiMedia Code : GRM051
Hydroxy-L-proline	: Corresponding HiMedia Code : RM052
L-Isoleucine	: Corresponding HiMedia Code : GRM053
L-Leucine	: Corresponding HiMedia Code : GRM054
L-Lysine monohydrochloride	: Corresponding HiMedia Code : GRM055
L-Methionine	: Corresponding HiMedia Code : GRM200
L-Phenylalanine	: Corresponding HiMedia Code : GRM060
L-Proline	: Corresponding HiMedia Codes : RM061
L-Serine	: Corresponding HiMedia Code : GRM063
L-Threonine	: Corresponding HiMedia Code : RM065
L-Tryptophan	: Corresponding HiMedia Code : GRM067
L-Tyrosine	: Corresponding HiMedia Code : RM069
L-Valine	: Corresponding HiMedia Codes : GRM070
Results	: All amino acids were tested and found to be passing as per standard specifications.

**Risk and Safety Information**

Storage Temperature(°C)	: Store below 30°C
-------------------------	--------------------

**Transport Information**

ADR/RID	: Not Dangerous Goods
---------	-----------------------

IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods