

---

## Potassium sulphate, Hi-LR™

**GRM404**

### Product Identifier

---

CAS No.	:	7778-80-5
EC No.	:	231-915-5
Molecular Formula	:	K <sub>2</sub> SO <sub>4</sub>
Molecular Weight	:	174.26
HS Code	:	3104 30 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Colourless or white to off-white crystals or powder or solid or granules
Solubility	:	33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	5.00 - 8.50
Chloride (Cl)	:	<= 0.05%
Iron (Fe)	:	<= 0.001%
Assay (Cation exchange titration)	:	min. 98.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TT5900000
WGK	:	1