

---

## Calcium propionate

**GRM4109**

### Product Identifier

---

CAS No.	:	4075-81-4
EC No.	:	223-795-8
Molecular Formula	:	C <sub>6</sub> H <sub>10</sub> CaO <sub>4</sub>
Molecular Weight	:	186.22
HS Code	:	2915 50 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Water (K.F.)	:	<= 9.50%
Assay (NT/EDTA T, on anhydrous basis)	:	97.00 - 100.50%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	UE8450000
WGK	:	2