
Sodium acetate trihydrate, Hi-LR™

GRM411

Product Identifier

CAS No.	:	6131-90-4
EC No.	:	204-823-8
Molecular Formula	:	C ₂ H ₃ O ₂ Na.3H ₂ O
Molecular Weight	:	136.08
HS Code	:	2918 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	7.50 - 9.20
Chloride (Cl)	:	<= 0.005%
Iron (Fe)	:	<= 0.002%
Sulfate (SO ₄)	:	<= 0.02%
Assay (NT)	:	99.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AJ4580000
WGK	:	1