
Sodium formate, Hi-AR™/ACS

GRM414

Product Identifier

CAS No.	:	141-53-7
EC No.	:	205-488-0
Molecular Formula	:	HCOONa
Molecular Weight	:	68.01
HS Code	:	2915 12 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white to off-white deliquescent crystals or powder or chunk
Solubility	:	100 mg soluble in 1 mL of water
Insoluble matter	:	<= 0.005%
Chloride (Cl)	:	<= 0.001%
Sulphate (SO ₄)	:	<= 0.001%
Calcium (Ca)	:	<= 0.005%
Iron (Fe)	:	<= 0.0005%
Heavy metals (as Pb)	:	<= 0.0005%
Assay (Iodometry)	:	99.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	LR0350000
WGK	:	1