
Sodium sulphite anhydrous, Hi-LRTM

GRM420

Product Identifier

CAS No.	:	7757-83-7
EC No.	:	231-821-4
Molecular Formula	:	Na ₂ SO ₃
Molecular Weight	:	126.04
HS Code	:	2832 20 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow crystals or powder or solid
Solubility	:	100 mg soluble in 1 mL of water
Chloride (Cl)	:	<= 0.05%
Iron (Fe)	:	<= 0.002%
Assay (Na ₂ S ₂ O ₃ Titration)	:	96.00 - 100.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WE2150000
WGK	:	1