
Isophthalic acid

GRM4318

Product Identifier

CAS No.	:	121-91-5
EC No.	:	204-506-4
Molecular Formula	:	C ₈ H ₆ O ₄
Molecular Weight	:	166.13
HS Code	:	2917 39 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder
Solubility	:	10 mg soluble in 10 mL of alcohol
FTIR	:	Matches with the standard pattern
Assay (NaOH Titration/GC)	:	99.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	NT2007000
WGK	:	2