
Ferrous sulphate, dried

GRM4376

Product Identifier

CAS No.	:	7720-78-7
Molecular Formula	:	FeSO ₄ .xH ₂ O
Molecular Weight	:	151.91 (anhydrous basis)
HS Code	:	2833 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow to buff or grey powder
Solubility	:	33.3 mg soluble in 1 mL of dilute sulphuric acid
Insoluble substances	:	<= 0.05%
Assay (KMnO ₄ /Ceric sulphate Titration)	:	86.00 - 89.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available