

---

## Silver bromide

**GRM4478**

---

### Product Identifier

CAS No.	:	7785-23-1
EC No.	:	232-076-8
Molecular Formula	:	AgBr
Molecular Weight	:	187.77
HS Code	:	2843 29 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	Yellow to yellow-green crystals or powder or chunks (darkens on exposure to light)
Solubility	:	Insoluble in water
Assay	:	min. 99.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3