

GRM4625 4-Chloroaniline**Product Number Packing**

GRM4625 : 500G

Product Information

Product Code : GRM4625
Product Name : 4-Chloroaniline
Molecular Formula : C₆H₆ClN
Molecular Weight : 127.57
CAS No. : 106-47-8
EC No. : 203-401-0
HS Code : 2921 42 00
EC Index No. : 612-137-00-9
Shelf Life : 4 years

Technical Specification

Appearance : White to brown crystals or powder or chips or chunks
Solubility : 100 mg soluble in 1 mL of alcohol
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 67 - 73°C
Assay (NT/GC) : 98.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H301-H311-H317-H331-H350-H410
Precautionary Statement(s) : P201-P261-P273-P280-P301+P310-P311
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 45-23/24/25-43-50/53
S-Phrase(s) : 53-45-60-61
WGK : 3
RTECS : BX0700000
Flash Point(°F) : 248°F
Flash Point(°C) : 120°C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 2018



Product Information

Class	:	6.1
Packaging Group	:	2
Marine Pollutant	:	Yes
ADR/RID	:	2018 6.1/PG 2
IMDG	:	2018 6.1/PG 2
IATA	:	2018 6.1/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein