

## Nickel ammonium sulphate hexahydrate, Hi-LR<sup>TM</sup>

GRM479


### Product Identifier

CAS No.	:	7785-20-8
EC No.	:	239-793-5
Molecular Formula	:	(NH <sub>4</sub> ) <sub>2</sub> Ni(SO <sub>4</sub> ) <sub>2</sub> ·6H <sub>2</sub> O
Molecular Weight	:	394.99
Synonym	:	Ammonium nickel sulphate hexahydrate, Hi-LR <sup>TM</sup>
HS Code	:	2833 24 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Green to blue-green crystals or powder or chunks
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (EDTA Titration)	:	min. 97.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302+H332- H317- H334- H341- H350i- H360D- H372- H410
Precautionary Statement(s)	:	P201- P260- P273- P280- P284- P304 + P340
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	WS6061000
WGK	:	3