
Sodium methoxide, 25-30% soln. in methanol

GRM4801


Product Identifier

CAS No.	:	124-41-4
Molecular Formula	:	CH ₃ ONa
Molecular Weight	:	54.02
HS Code	:	3822 00 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow liquid
Solubility	:	1 mL miscible in 1 mL of ethanol
Assay (GC)	:	25.00 - 30.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H225- H251- H301+H311+H331- H314- H370
Precautionary Statement(s)	:	P210.- P235 + P410- P260- P280- P301 + P310- P305+P351+P338
UN No.	:	1289
Class	:	3
Packing Group	:	II
RTECS	:	Not available
WGK	:	1