
Taurine, Hi-LRTM

GRM4806

Product Identifier

CAS No.	:	107-35-7
EC No.	:	203-483-8
Molecular Formula	:	C ₂ H ₇ NO ₃ S
Molecular Weight	:	125.15
HS Code	:	2921 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder.
Solubility	:	33.3 mg soluble in 1 mL of water
pH (4% in water at 25°C)	:	4.50 - 6.00
FTIR	:	Matches with the standard pattern
Assay (NaOH Titration)	:	99.00 - 103.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WX0175000
WGK	:	2