

**GRM4815 Tetramethylammonium bromide****Product Number Packing**

|         |        |
|---------|--------|
| GRM4815 | : 100G |
| GRM4815 | : 500G |


**Product Information**

|                   |                                      |
|-------------------|--------------------------------------|
| Product Code      | : GRM4815                            |
| Product Name      | : Tetramethylammonium bromide        |
| Molecular Formula | : C <sub>4</sub> H <sub>12</sub> BrN |
| Molecular Weight  | : 154.05                             |
| CAS No.           | : 64-20-0                            |
| EC No.            | : 200-581-2                          |
| HS Code           | : 2923 90 00                         |
| Shelf Life        | : 4 years                            |


**Technical Specification**

|                 |  |
|-----------------|--|
| Appearance      | : White to light grey crystals or powder |
| Solubility      | : 33.3 mg soluble in 1 mL of water       |
| FTIR (KBr disc) | : Matches with the standard pattern      |
| Assay (AT)      | : 98.00 - 102.00%                        |

**GHS Safety Information**

|                            |   |
|----------------------------|---|
| Hazard Statement(s)        | : H300-H315-H319-H335   |
| Precautionary Statement(s) | : P261-P264-P301+P310-P305+P351+P338  |
| Signal Word                | : Danger  |
| Hazard Pictogram(s)        |  |

**Risk and Safety Information**

|                         |   |
|-------------------------|---|
| R-Phrase(s)             | : 25-36/37/38   |
| S-Phrase(s)             | : 26-28-36-45   |
| WGK                     | : 3   |
| RTECS                   | : BS7600000   |
| Storage Temperature(°C) | : Store below 30°C  |
| Hazard Symbol(s)        |  |

**Transport Information**

|                  |                 |
|------------------|-----------------|
| UN No.           | : 2811          |
| Class            | : 6.1           |
| Packaging Group  | : 2             |
| Marine Pollutant | : No            |
| ADR/RID          | : 2811 6.1/PG 2 |



## Product Information

IMDG : 2811 6.1/PG 2  
IATA : 2811 6.1/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein