
Azure B, Hi-CertTM

GRM486

Product Identifier

CAS No.	:	531-55-5
EC No.	:	208-511-2
C.I. No.	:	52010
Molecular Formula	:	C ₁₅ H ₁₆ N ₃ ClS
Molecular Weight	:	305.83
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Green to brown to black crystals or powder
Solubility	:	20 ppm in water yields clear solution
Absorption maxima	:	640 - 655 nm (20 ppm in water)
Loss on drying (at 110°C, 2 hr)	:	<= 5.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	SO5550000
WGK	:	3