

## 4-Chlorobenzoic acid

**GRM504**


### Product Identifier

CAS No.	:	74-11-3
EC No.	:	200-805-9
Molecular Formula	:	C <sub>7</sub> H <sub>5</sub> ClO <sub>2</sub>
Molecular Weight	:	156.57
Synonym	:	p-Chlorobenzoic acid
HS Code	:	2916 39 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to yellow crystals or powder
Solubility	:	10 mg soluble in 10 mL of methanol
FTIR	:	Matches with the standard pattern
Melting range	:	235 - 245°C
Assay (GC/NaOH Titration)	:	98.00 - 101.50%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DG4976010
WGK	:	2