

Paraldehyde, Hi-LR™

GRM5376


Product Identifier

CAS No.	:	123-63-7
EC No.	:	204-639-8
Molecular Formula	:	C ₆ H ₁₂ O ₃
Molecular Weight	:	132.16
HS Code	:	2912 50 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow liquid
Solubility	:	1 mL miscible with 1 mL of alcohol
FTIR	:	Matches with the standard pattern.
Density (at 25°C)	:	0.989 - 0.999 g/mL
Assay (GC)	:	min. 97.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H226
Precautionary Statement(s)	:	P210- P233- P240- P241- P242- P243
UN No.	:	1264
Class	:	3
Packing Group	:	III
RTECS	:	YK0525000
WGK	:	1