

## Tetrabutylammonium bromide, Hi-AR®

**GRM1162**


### Product Identifier

CAS No.	:	1643-19-2
EC No.	:	216-699-2
Molecular Formula	:	C <sub>16</sub> H <sub>36</sub> BrN
Molecular Weight	:	322.37
HS Code	:	2923 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to slightly cream hygroscopic crystals or powder or granules
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Assay (AT)	:	99.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H361fd- H412
Precautionary Statement(s)	:	P201- P273- P280- P308 + P313- P337 + P313- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3