

p-Hydroxybenzoic acid, Hi-LRTM

GRM557


Product Identifier

CAS No.	:	99-96-7
EC No.	:	202-804-9
Molecular Formula	:	C ₇ H ₆ O ₃
Molecular Weight	:	138.12
Synonym	:	4-Hydroxybenzoic acid
HS Code	:	2918 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to cream crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of acetone
FTIR	:	Matches with the standard pattern
Melting range	:	212°C - 217°C
Assay (HPLC/GC/NaOH Titration)	:	98.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H318- H335
Precautionary Statement(s)	:	P280- P305+P351+P338+P310
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH1925000
WGK	:	1